# Handleiding secMAC

Video Laryngoscoop

CE 0197

Model: E.An-IILP & E.An-IISP

**Tianjin Medan Medical Corp.**

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8.6 Color:  $\geq 5000$ k

8.7 Anti-fog function

8.8 Imaging display, taking photos, taking video

8.9 Imaging storage: volume 4G/8G

8.10 IP44

## 9. Specifications of external power

INPUT:AC100-240V~ 50/60HZ,OUTPUT:DC5V 1000mA

## 10. Environments

10.1 Working environments:

- Temperature:  $10^\circ\text{C} \sim 40^\circ\text{C}$
- Humidity :  $30\% \sim 75\%$
- Atmospheric pressure: 700hpa  $\sim$  1060hpa
- Adaptation power: DC 4.2V (external adapter voltage: DC5V,1000 mA)

10.2 Transportation, storage environments

- Temperature:  $-40^\circ\text{C} \sim +55^\circ\text{C}$
- Humidity:  $\leq 93\%$
- Atmospheric pressure: 500hpa  $\sim$  1060hpa

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## 12. Electromagnetic Compatibility

The Laryngoscope is suitable for use in an electromagnetic environment as described in this chapter.

### CAUTION

The instrument has been tested and complies with the limits of IEC 60601-1-2. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices which can be determined by turning the equipment off and on, the user

is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device.
- Increase the separation between the equipment.
- Consult the manufacturer for help.

Table 12.1

Guidance and manufacturer's statement - Electromagnetic emission		
Laryngoscope is intended to be used in the following specified electromagnetic environment, the customer or the user should ensure that it is used in such electromagnetic environment:		
Emission Test	Compliance	Electromagnetic Environment - Guide
RF Emissions IEC/CISPR11:2010	1 item	The laryngoscope uses RF energy only for its internal function. Therefore, its RF emissions are very low are not likely to cause any interference in nearby electronic equipment.
RF Emissions IEC/CISPR11:2010	Type B	The laryngoscopes are suitable for use in all facilities that are not directly connected to non-domestic and public residential low-voltage supply networks.
Harmonic emission IEC 61000-3-2-2009	No application	
Voltage fluctuations/ flicker emissions IEC61000-3-3:2005	No application	

Table 12.2

Guidance and manufacturer's statement - Electromagnetic Immunity			
Laryngoscope is intended to be used in the following specified electromagnetic environment, the customer or the user should ensure that it is used in such electromagnetic environment:			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment: Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	+/- 6 kV contact +/- 8 kV air	+/- 6 kV contact +/- 8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical Fast transient/burst IEC 61000-4-4	+/- 2 kV for power supply line	+/- 2 kV for power supply line	The network power supply should have a typical commercial or hospital environment quality
Surge IEC 61000-4-5	+/- 1 kV differential mode +/- 2 kV common mode	+/- 1 kV differential mode +/- 2 kV common mode	The network power supply should have a typical commercial or hospital environment quality
Voltage dips, short interruptions and voltage variations on power supply	<5% U T (>95% dip in U T) for 0.5 cycle 40% U T (60% dip in U T) for 5 cycles 70% U T (>30% dip in U T) for 25 cycles <5% U T (>95% dip in U T) for 5	<5% U T (>95% dip in U T) for 0.5 cycle 40% U T (60% dip in U T) for 5 cycles 70% U T (>30% dip in U T) for 25 cycles <5% U T (>95% dip in U T) for 5	The network power supply should have a typical commercial or hospital environment quality. If the user of the video laryngoscope requires continuous operation during a power interruption, it is recommended that the video laryngoscope be powered by

input lines IEC 61000-4-11	seconds	seconds	uninterruptible power supply or battery
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	The Power Frequency Magnetic Field should be of typical commercial or hospital environment quality
NOTE: U T is the AC mains voltage prior to application of the test level.			

Table 12.3

Guidance and manufacturer's statement - Electromagnetic immunity			
Laryngoscope is intended to be used in the following specified electromagnetic environment, the customer or the user should ensure that it is used in such electromagnetic environment:			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment: Guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 V(Valid values) 150 kHz to 80 MHz	3 V(Valid values)	Portable and mobile RF communications equipment should be used no closer to any part of the laryngoscope, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P} \quad 80 \text{ MHz} \sim 800 \text{ MHz}$ $d = 2.3\sqrt{P} \quad 800 \text{ MHz} \sim 2.5 \text{ GHz}$ In formula: Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and D is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
RF Irradiated IEC 61000-4-3	3 V/m 80 MHz~2.5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from surfaces, objects, and people.

- a) Fixed-range transmitters, such as wireless (cellular / cordless) telephones and terrestrial mobile radio base stations, amateur radio, AM and FM radio and television broadcasts, are not reasonably predictable in theory. In order to evaluate the electromagnetic environment of fixed RF transmitters, the investigation of electromagnetic sites should be considered. If the field strength of the site where the video laryngoscope is measured is higher than the applicable RF compliance level, the video laryngoscope should be observed to verify that it is operating normally. If abnormal performance is observed, then additional action may be necessary, such as reorienting the direction or position of the video laryngoscope.
- b) The field strength should be less than 3 V / m in the entire frequency range from 150 kHz to 80 MHz.

Table 12.4

Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the laryngoscope			
<p>The laryngoscope is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The user of the laryngoscope can help to prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the laryngoscope as recommended below, according to the maximum output power of the communications equipment.</p>			
Rated Maximum Output Power of Transmitter (W)	Separation Distance According to Frequency of Transmitter (m)		
	150 kHz~80 MHz $d = 1.2\sqrt{P}$	80 MHz~800 MHz $d = 1.2\sqrt{P}$	800 MHz~2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
<p>For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p> <p>NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by the absorption and reflection from structures, objects, and people.</p>			